

Instrumentation for Linac-based X-Ray FELs

Henrik Loos

12th Beam Instrumentation Workshop
May 1-4, 2006

Outline

- XFEL introduction
- LCLS overview
- Electron beam diagnostics
 - Transverse Beam Properties
 - Longitudinal Beam Properties
- Photon beam diagnostics

X-Ray FEL Features

- ~1Å photon wavelength or ~10keV photon energy
- Uses SASE principle to amplify and saturate spontaneous radiation in ~100m of undulator
- Requires
 - Multi GeV beam energy
 - kA peak beam current
 - Micron beam emittance to match photon beam phase space


$$\lambda_R = \frac{\lambda_U}{2\gamma^2} (1 + K_{rms}^2)$$

$$\varepsilon \approx \frac{\lambda}{4\pi}$$


X-Ray FEL Parameters

Electron Beam		LCLS	XFEL	SCSS
Energy	GeV	4.3-13.6	10-20	6.1
Peak Current	kA	3.4	5	3
Bunch Charge	nC	0.2-1	1	1
Norm. Slice Emittance	μm	1.2	1.4	0.85
Bunch Length	fs	70	80	80
Slice Energy Spread	MeV	1.4	2.5	0.25
Photon Beam		LCLS	XFEL	SCSS
Saturation Length	m	60-100	40-170	80
Photon Energy	keV	0.8-8	0.2-12.4	12
Peak Power	GW	4-8	22-135	3

Linac Coherent Light Source



LCLS Accelerator Layout



Courtesy P. Emma

3 May 2006

BIW 2006

Henrik Loos

loos@slac.stanford.edu


LCLS Diagnostics Tasks

- Charge
 - Toroids (Gun, Inj, BC, Und)
 - Faraday cups (Gun & Inj)
- Trajectory & energy
 - Stripline BPMs (Gun, Inj, Linac)
 - Cavity BPMs (Und)
 - Profile monitors (Inj), compare position with alignment laser
- Transverse emittance & energy spread
 - Wire scanners
 - YAG screen (Gun, Inj)
 - OTR screens (Inj, Linac)
- Bunch length
 - Cherenkov radiators + streak camera (Gun)
 - Transverse cavity + OTR (Inj, Linac)
 - Coherent radiation power (BC)
- Slice measurements
 - Horizontal emittance
 - T-cavity + quad + OTR
 - Vertical Emittance
 - OTR in dispersive beam line + quad
 - Energy spread
 - T-cavity + OTR in dispersive beam line


Diagnostics Requirements

Parameter	Method	Unit	Resolution
Current	Toroid, FC	%	2
Position	Stripline BPM	μm	5 - 20
	Cavity BPM	μm	1
Beam Size	Wire Scanner	μm	5
	YAG	μm	15 – 30
	OTR	μm	5 – 30
Bunch Length	Streak Camera	fs	300
	Transverse Cavity	Slices	10
	BLM	%	5

LCLS Injector Diagnostics



LCLS Linac Diagnostics



SLAC linac tunnel

research yard

Beam Profile Monitors (YAG & OTR)

■ YAG requirements

- Use 100μm thick crystals to meet resolution
- GTF measurements show feasibility

■ OTR requirements

- Optimize yield to enable beam profile measurement at 0.2nC


OTR yield for 100mrad angular acceptance

Energy (MeV)	QE (%), 450-650 nm	QE (%), 400-750 nm
135	0.44	0.75
4300	0.98	1.68
13500	1.17	1.99


- Provide sufficient depth of field for imaging of 45° foil
 - Simulation shows 1mm DOF for f/# of 5 within 20μm resolution
- Match direction of reflection with axis of dispersion or T-CAV deflection
- Foil is aluminum to optimize TR yield and 1μm thick to minimize radiation

Optics Layout

- Used for all standard YAG/OTR screens
- Telecentric lens
 - 55mm focal length
 - >100 line pairs/mm
 - Magnification up to 1:1
- Stack of 2 insertable neutral density filters
- Beam splitter and reticule for in situ calibration
- Megapixel CCD with 12bit and 4.6 μ m pixel
- Radiation shielding required in main linac tunnel




OTR/YAG Optics Design



Courtesy V. Srinivasan

OTR Imager with Tilted Geometry


- Need wide field of view in focus for measurements in spectrometer beam line
- Tilt OTR screen and CCD by 5 degrees in 1:1 imaging
- 10um resolution



B.X. Wang et al. PAC05

Simulation of OTR Beam Size Measurement

- Simulation of CCD image
 - Include 0.5% TR yield, photon shot noise, and typical CCD parameters for quantum efficiency, read out noise, pixel size, digitizer gain
- Calculation of beam size
 - Generate beam profile with 10σ bounding box
 - Compare rms width of profile with original Gaussian beam size
- Simulation agrees well with OTR measurement at GTF
- Error of 5% in beam size for beam of 0.1nC, 260 μ m at 10 μ m resolution



Longitudinal Diagnostics

- Gun region
 - Cherenkov radiator & streak camera
- Bunch length and slice emittance
 - Transverse cavity
- Longitudinal feedback loop
 - Integrated power from coherent radiation

Cherenkov Radiators

- Located in gun region for temporal diagnostics of 6 MeV beam from gun
- Convert electron beam time structure into light pulse for streak camera measurement
- Cherenkov light suitable at low beam energies
- Design requirements
 - Match time resolution of radiator to streak camera (Hamamatsu FESCA-200, $< 300\text{fs}$)
 - Generate and transport a sufficient # of photons for 200pC beam to streak camera in laser room (10m away)


Cherenkov Radiator Design

■ Fused silica

- $n = 1.458$, $\theta_{CR} = 46.7^\circ$
- Total internal reflection
- Frosting of back surface
- $N_\phi = 7.5/e/mm/50nm$
@400nm

■ Temporal and spatial resolution


- Thickness of $100\mu m$
- $\Delta t = 375fs$
- $\Delta x = 190\mu m$



Courtesy D. Dowell


Optical Transport Layout

- 1:1 relay imaging from radiator to streak camera
- Assume 1% efficiency from frosting to scatter into 100mrad
- 6% acceptance through tube for source of 5mm x 100mrad
- $1.5 \cdot 10^5$ photons on slit of streak camera for 200 pC



Courtesy D. Dowell


Transverse Cavity




- Translates longitudinal into transverse beam profile when operating at RF zero crossing
- Parasitic operation with kicker and off-axis screen
- Single shot absolute bunch length measurement
- Temporal resolution limited by unstreaked spot size

Transverse Cavity Measurement at TTF

Beam without and with BC 3 (second bunch compressor)



1 picosecond




1 picosecond

Scans at high power ~16MW

Courtesy J. Frisch

TCAV in LCLS after BC2

- Short 70fs bunch length requires full RF power for cavity
- Parasitic measurement with beam optics optimized for SASE
- Resolution 20fs sufficient for length measurement





Bunch Length Monitor

- Relative bunch length measurement used for longitudinal feedback
- Non-intercepting, calibrated with interceptive TCAV measurement
- Based on integrated power from coherent radiation source (C*R)

$$W = N_e^2 \int d\omega \frac{dW_1(\omega)}{d\omega} f(\omega), \quad f(\omega) = \left| \int n(t) e^{i\omega t} dt \right|^2$$

- Single electron radiation spectrum $W_1(\omega)$ depends on radiation source
- Bunch length determined by long wavelengths $\lambda \gg 2\pi\sigma_{\text{rms}}$
- BC1: 1cm – 1mm
- BC2: 1mm - .1mm




Radiation Sources

- Wide range of bunch lengths from 25um to 300um
- Diode detectors work well below 300GHz
- Pyroelectric detectors work well above 300GHz
- Long bunches
 - Couple radiation from ceramic gap in beam pipe into waveguides with different diode detectors
- Short bunches
 - Extract coherent radiation from bend magnet with hole mirror and send to a pyroelectric detector



CER Detector Layout

- Edge rad. dominates over synchrotron and diffraction
- Near field calculation necessary for radiation spectrum at detector




Bunch Length Sensitivity of Detector Signal

- Detection efficiency includes diffraction, vacuum window, water absorption, pyroelectric detector response, and bunch form factor.
- Introduce high and low pass filters at 10cm^{-1} and 20cm^{-1} .



Gap Radiation Detector

- Expect 2uJ radiation energy from 2cm gap for 1nC, 200um bunch
(Calculation J. Wu)
- Energy density of 1.6nJ/mm^2
- Diode sensitivity $\sim 0.1\text{pJ/mm}^2$
- Disperse pulse in 20cm waveguide to keep diodes in linear range
- Diodes paired to reduce dependence on beam position




Courtesy S. Smith

Photon Beam Diagnostics

- Measure spontaneous radiation for undulator commissioning
- Measure FEL photon beam for SASE commissioning
- Nondestructive measurements of beam properties for user operation

LCLS FEE Schematic



Wide Field of View Direct Imager




CAMERA
Link


Single shot measurement of
 $f(x,y)$, x , y , u





Camera



Scintillators




Photoelectrons
generated by 0.01%
FEL



Courtesy R. Bionta


Indirect Imager

Single shot measurement of $f(x,y)$, x , y , u
Multi shot measurement of λ




Material Pairs	N	d	Γ	Primary Energy
B ₄ C/SiC	150	60 Å	0.7	8.261 keV
B ₄ C/SiC	35	100 Å	0.75	8.261 keV
Be/SiC	40	110 Å	0.75	0.8261 keV
B ₄ C/SiC	750	20 Å	0.55	24.78 keV
B ₄ C/SiC	450	20 Å	0.55	24.78 keV
B ₄ C/SiC	150	70 Å	0.45	24.78 keV
B ₄ C/SiC	50	100 Å	0.754	2.478 keV

B4C/SiC Test Multilayers Fabricated




Angle selects energy and attenuation



Courtesy R. Bionta


Total Energy Calorimeter

Single shot measurement of
 $f(x,y)$, x , y , u




Thermal diffusion calculations
performed

$t = 300 \mu\text{s}$ $t = 100 \mu\text{s}$




$t = 0$




Courtesy R. Bionta

Ion Chamber

Single shot, *non destructive*,
measurement of
 x', y', x, y, u



Segmented
cathodes for
position
measurement



Imaging of optical
emission for
position
measurement

Courtesy R. Bionta

Summary

- Electron beam diagnostics based on proven methods
- Photon beam diagnostics needs development of new techniques which are difficult to test due to the lack of a photon source comparable to an X-FEL

Acknowledgements

- Thanks to many colleges from the LCLS collaboration