

Using IDL and Python with EPICS

Mark Rivers, University of Chicago

Argonne National Laboratory



*A U.S. Department of Energy
Office of Science Laboratory
Operated by The University of Chicago*



Outline

- **Quick Overview of IDL**
- **ezca library**
- **Calling ezca from IDL**
- **IDL CA API**
- **IDL EPICS class libraries**
- **IDL applications**
- **Overview of Python**
- **Python class libraries**
- **Python applications**
- **Using EPICS from Visual Basic**

Overview of IDL

- A high-level interpreted *programming language* with vector and array primitives - sort of a cross between BASIC and APL
- Modern programming language
 - Flow control
 - Data structures
 - Objects
- All operators and most functions work on scalar, vector or array data of any data type.
- Data *visualization tool*, advanced built-in graphics
 - 2-D plots
 - Contour plots
 - Surface plots
 - Shaded surfaces
 - Gray scale/color images
 - Isosurfaces
 - Volume rendering
- Multi-platform support
 - Unix: Sun, Hewlett Packard, Silicon Graphics, IBM
 - Linux
 - Microsoft Windows
 - Mac Darwin
- List price: ~\$3,000 on workstations, ~\$1,500 on PC, Mac

Overview of IDL

- **Can call external C or other code**
- **Very fast for array operations, as fast as compiled languages**
- **GUI builder**
- **Multi-threaded**
- **Good vendor support**
- **IDL GUI applications can be run at no cost (IDL 6.0 and above)**
 - **Must have license to use IDL command line**

Overview of IDL

Data Structures

- **A variable in IDL has both a structure and a data type associated with it. Both of these are dynamic, i.e. they can be changed via an assignment statement at any time.**
- **Data types**
 - Byte (b bit, unsigned)
 - Integer (16 bit, signed)
 - Long (32 bit, signed)
 - Float (32 bit floating point)
 - Double (64 bit floating point)
 - Complex (pair of 32 bit floats)
 - Double complex (pair of 64 bit floats)
 - String (0 to 64k characters)
- **Data Structures**
 - Scalar
 - Vector
 - Array - up to 7 dimensions
 - Structure - composed of other elements, like C
- **The sizes of arrays are limited only by the available virtual memory.**

Overview of IDL

Assignment Statements

`A = B + 1`

- **A has the same structure as B, with a data type equal to that of the most precise operand in the expression on the right hand side. In this case it could be any type except string.**
- **If B is a vector or array then 1 is added to each element.**

`A = 0` ; A is a 16 bit integer

`A = A * 0.5` ; A is now a 32 bit float

`B = A(*,3)` ; B is equal to the 4th row of A

`A(*,3) = 0` ; Set all elements in 4th row of A equal to 0

Syntax

- **Examples:**

`image = fltarr(512, 512)` ; zero filled array

`b = image(0:127, 0:127)` ; b is 128x128 array

`image(*,100) = findgen(512)` ; replace row 100

`plot, image(*,120)` ; plot row 121

; Display the power spectrum as an image

`tvscrl, alog(abs(fft(image, 1)))`


IDL Examples

```
IDL> a = sin(findgen(100)/99. * 2 * !pi)
```

```
IDL> help, a
```

```
A                FLOAT      = Array[100]
```


```
IDL> plot, a
```



IDL Examples

```
IDL> a = shift(alog(abs(fft(dist(256),1))),128,128)
```

```
IDL> isurface, a
```



ezca and EzcaScan

- **EPICS extensions for “Easy Channel Access”**

- Don't need to handle chids, just use PV name strings – hash table
- Synchronous APIs - applications don't have to handle callbacks

- **Ezca (partial list)**

```
- epicsShareFunc int epicsShareAPI ezcaGet(char *pvname, char ezcatype,  
-     int nelem, void *data_buff);  
- epicsShareFunc int epicsShareAPI ezcaPut(char *pvname, char ezcatype,  
-     int nelem, void *data_buff);  
- epicsShareFunc int epicsShareAPI ezcaPutOldCa(char *pvname, char ezcatype,  
-     int nelem, void *data_buff);  
- epicsShareFunc int epicsShareAPI ezcaNewMonitorValue(char *pvname,  
-     char ezcatype); /* returns TRUE/FALSE */  
- epicsShareFunc int epicsShareAPI ezcaSetTimeout(float sec);  
- epicsShareFunc float epicsShareAPI ezcaGetTimeout(void);  
- epicsShareFunc int epicsShareAPI ezcaSetRetryCount(int retry);  
- epicsShareFunc int epicsShareAPI ezcaGetRetryCount(void);  
- epicsShareFunc int epicsShareAPI ezcaPvToChid(char *pvname, chid **cid);  
- epicsShareFunc int epicsShareAPI ezcaSetMonitor(char *pvname, char ezcatype);  
- epicsShareFunc int epicsShareAPI ezcaClearMonitor(char *pvname, char ezcatype);  
- epicsShareFunc int epicsShareAPI ezcaStartGroup(void);  
- epicsShareFunc int epicsShareAPI ezcaEndGroup(void);  
- epicsShareFunc int epicsShareAPI ezcaGetControlLimits(char *pvname,  
-     double *low, double *high);  
- epicsShareFunc int epicsShareAPI ezcaGetGraphicLimits(char *pvname,  
-     double *low, double *high);  
- epicsShareFunc int epicsShareAPI ezcaGetNelem(char *pvname, int *nelem);  
- epicsShareFunc int epicsShareAPI ezcaGetPrecision(char *pvname,  
-     short *precision);  
- epicsShareFunc int epicsShareAPI ezcaGetStatus(char *pvname,  
-     TS_STAMP *timestamp, short *status, short *severity);  
- epicsShareFunc int epicsShareAPI ezcaGetUnits(char *pvname,  
-     char *units); /* units must be at least EZCA_UNITS_SIZE large */
```

ezca and EzcaScan

- **EzcaScan (partial list)**

- `epicsShareFunc int epicsShareAPI Ezca_getArray(int noName, char **pvName, int type, int noda, void *value);`
- `epicsShareFunc int epicsShareAPI Ezca_getArrayEvent(int noName, char **pvName, int type, int noda, void *value);`
- `epicsShareFunc int epicsShareAPI Ezca_putArray(int noName, char **pvName, int type , int noda, void *value);`
- `epicsShareFunc int epicsShareAPI Ezca_putArrayEvent(int noName, char **pvName, int type, int noda, void *value);`

ezca and IDL

- **IDL can call “shareable libraries”, e.g. .so files on Unix, .dll files on Windows**
- **The argument passing convention is fixed, it is not compatible with ezca.dll directly**
- **Need a thin glue layer between IDL and ezca/EzcaScan**
- **ezcaIDL is the glue layer. Mostly just changes calling conventions. Provides a few functions that ezca and EzcaScan do not. Use the ezcaPvToChid() function.**
 - ezcaIDLGetEnumStrings()
 - ezcaIDLGetCountAndType()

IDL Channel Access API

Routines which return information about process variables

```
Status = caGet(pvname, value, /string, maximum=max)
Status = caGetControlLimits(pvname, low, high)
Status = caGetGraphicLimits(pvname, low, high)
Status = caGetPrecision(pvname, precision)
Status = caGetStatus(pvname, timestamp, status, severity)
Status = caGetUnits(pvname, units)
Status = caGetEnumStrings(pvname, strings)
Status = caGetCountAndType(pvname, count, type)
```

Routines which write new values to process variables

```
Status = caPut(pvname, value, wait=wait)
```

Routines which control channel access timeouts

```
Timeout = caGetTimeout()
caSetTimeout, timeout
RetryCount = caGetRetryCount()
caSetRetryCount, retrycount
```

IDL Channel Access API

Routines which control synchronous groups

`caStartGroup`

`stat = caEndGroup(status)`

Routines which control channel access monitors

`Status = caSetMonitor(pvname)`

`Status = caClearMonitor(pvname)`

`State = caCheckMonitor(pvname)`

Routines which control debugging and error messages

`caDebug, state`

`caTrace, state`

`caError, err_string, /ON, /OFF, /PRINT, prefix=prefix`



Documentation at

<http://cars.uchicago.edu/software/idl/ezcaIDLGuide.html>

<http://cars.uchicago.edu/software/idl/ezcaIDLRef.html>

IDL EPICS Examples

```
IDL> status = caget('13LAB:m1.VAL', position)
IDL> help, status, position
STATUS          LONG          =          0
POSITION        DOUBLE        =      517.19305
IDL> status = caget('13LAB:quadEM:mca1', spectrum)
IDL> plot, spectrum
IDL> help, status, spectrum
STATUS          LONG          =          0
SPECTRUM        LONG          = Array[2048]
IDL> plot, spectrum[0:500]
IDL> fft_data = alom(abs(fft(spectrum,1)))
IDL> plot, fft_data, xrange=[0,1023]
```



IDL EPICS Examples

Move a motor

```
DL> status = caput('13LAB:m8.VAL', 10000)
IDL> status = caget('13LAB:m8.RBV', pos)
IDL> print, pos
      215.52734
IDL> status = caget('13LAB:m8.RBV', pos)
IDL> print, pos
      835.64453
IDL> status = caget('13LAB:m8.RBV', pos)
IDL> print, pos
     1795.6055
```

Ezca timeout values are important!

```
IDL> print, cagettimeout()
      0.0200000
IDL> t0=sysime(1)&for i=1,100 do t=caget('13LAB:m1', v)&print, sysime(1)-t0
      2.9898720
IDL> casetimeout, .001
IDL> t0=sysime(1)&for i=1,100 do t=caget('13LAB:m1', v)&print, sysime(1)-t0
      0.21649790
```

IDL EPICS Examples

Using monitors

Monitored channels read the cached values on caget()

Can check whether a monitor has happened (a Channel Access value callback)

```
IDL> status = caSetMonitor('13LAB:m8.DMOV')
IDL> state = caCheckMonitor('13LAB:m8.DMOV')
IDL> help, state
STATE          LONG          =          1
IDL> status = caget('13LAB:m8.DMOV', done)
IDL> help, done
DONE           INT           =          1
IDL> state = caCheckMonitor('13LAB:m8.DMOV')
IDL> help, state
STATE          LONG          =          0
IDL> status = caput('13LAB:m8.VAL', 0)
IDL> state = caCheckMonitor('13LAB:m8.DMOV')
IDL> help, state
STATE          LONG          =          1
IDL> status = caget('13LAB:m8.DMOV', done)
IDL> help, state
STATE          LONG          =          1
IDL> help, done
DONE           INT           =          0
```

Monitors are useful for seeing that a PV changed state, even if its value is the same because one “missed” the transition. For example, PV goes 0->1->0. IDL polling might miss the one state, but checking a monitor would let one know that it happened.

IDL EPICS Class Libraries

- **IDL object classes that hide the underlying EPICS process variables**
- **IDL objects treat all data as private, only accessible through methods.**
- **Provide an object-oriented interface to common beamline objects (motors, scalers, mcas, scans)**
 - `epics_motor`
 - `epics_scaler`
 - `epics_mca` (inherits device-independent mca class)
 - `epics_med` (multi-element detector)
 - `epics_sscan`
- **Example of `epics_motor`**

```
IDL> motor = obj_new('EPICS_MOTOR', '13LAB:m8')
```


```
IDL> motor->move, 1000.      ; Move to absolute position 10.
```

```
IDL> motor->wait              ; Wait for it to get there
```

IDL EPICS Class Libraries

- **Example of epics_mca**

```
IDL> mca = obj_new('epics_mca', '13LAB:aim_adc1')
IDL> mca->erase
IDL> mca->acquire_on
IDL> data = mca->get_data()
IDL> plot, data
```



- **Example of epics_scaler**

```
IDL> scaler = obj_new('epics_scaler', '13LAB:scaler1')
IDL> scaler->start, 10.      ; Count for 10 seconds
IDL> scaler->wait            ; Wait for it to get done
IDL> counts = scaler->read(); Read the counts on each channel
IDL> print, counts
100000000      0      0      0      0      0      0      0
```

IDL EPICS Class Libraries


- **epics_sscan**
- **Designed to do the following:**
 - Provide an object-oriented interface to standard EPICS scans, enabling user written software to easily access scan header information and data.
 - Provide an easy way to read MDA files written by the saveData function in synApps.
 - Provide an easy way to get scan data into the IDL iTools system. iTools provide powerful interfaces for visualizing data, zooming in, adding annotation, and producing publication quality plots.
 - Provide a way to convert binary scan files (e.g. MDA) into ASCII
 - Does not currently communicate with the IOC for real-time data, but this is planned for the future

IDL EPICS Class Libraries

Example: Simple 1D epics_sscan

```
IDL> s = read_mda('13IDC_0027.mda') ; Read the data
```


```
IDL> s->display ; Display the first detector
```



IDL EPICS Class Libraries

Example: 2-D epics_sscan

```
IDL> s=read_mda('2idd_0087.mda') ; Read the 2-D dataset  
IDL> s->display, /all, /grid ; Display all of the  
images in a grid
```




IDL EPICS Class Libraries

Example: 2-D epics_sscan

; Plot a profile of column 20 (X=20) in detector 15.

```
IDL> s->display, detector=15, xrange=20
```



IDL EPICS Class Libraries

Documentation: Reference manual for each class library

EPICS_MCA Class

This page was created by the IDL library routine `mk_html_help`. For more information on this routine, refer to the IDL Online Help Navigator or type:

?mk_html_help at the IDL command line prompt.

Last modified: Sat Jul 14 10:16:05 2001.

List of Routines

- [EPICS_MCA::ACQUIRE OFF](#)
- [EPICS_MCA::ACQUIRE ON](#)
- [EPICS_MCA::ACQUIRE WAIT](#)
- [EPICS_MCA::ADD ROI](#)
- [EPICS_MCA::DEL ROI](#)
- [EPICS_MCA::ERASE](#)
- [EPICS_MCA::GET ACQUIRE STATUS](#)
- [EPICS_MCA::GET CALIBRATION](#)
- [EPICS_MCA::GET DATA](#)
- [EPICS_MCA::GET ELAPSED](#)
- [EPICS_MCA::GET PRESETS](#)
- [EPICS_MCA::GET ROIS](#)
- [EPICS_MCA::GET ROI COUNTS](#)
- [EPICS_MCA::GET SEQUENCE](#)
- [EPICS_MCA::INIT](#)
- [EPICS_MCA::SET CALIBRATION](#)
- [EPICS_MCA::SET DATA](#)
- [EPICS_MCA::SET PRESETS](#)
- [EPICS_MCA::SET ROIS](#)
- [EPICS_MCA::SET SEQUENCE](#)
- [EPICS_MCA::SPECTRA SCAN](#)
- [EPICS_MCA::WRITE FILE](#)
- [EPICS_MCA_DEFINE](#)
- [RELEASE NOTES](#)

Documentation at:

<http://cars9.uchicago.edu/software/idl/>

EPICS_MCA::GET_DATA

[\[Previous Routine\]](#) [\[Next Routine\]](#) [\[List of Routines\]](#)

NAME:

EPICS_MCA::GET_DATA

PURPOSE:

This function returns the data from the MCA.

CATEGORY:

EPICS device class library.

CALLING SEQUENCE:

Result = epics_mca->GET_DATA()

KEYWORD PARAMETERS:

CHECK_NEW:

A flag which indicates that this routine should only return the data if it has changed.

OPTIONAL OUTPUTS:

NEW_FLAG:

If CHECK_FLAG is set, then NEW_FLAG will be 1 if the function is returning new data, 0 if the function is not returning new data. If CHECK_FLAG is set and NEW_FLAG is 0 then the function returns -1.

PROCEDURE:

This function reads the data from the hardware using the EPICS MCA record, and then invokes MCA::GET_DATA

ADDITIONAL INFORMATION:

See [MCA::GET_DATA\(\)](#).



MODIFICATION HISTORY:

Written by: Mark Rivers, October 1, 1997
Nov. 14, 1997 Mark Rivers. Changed routine to eliminate setting rec.READ back to 0, since record support does this automatically and it was causing record to process again.
19-Sep-1998 MLR Added /WAIT to caput, since default is not to wait for callback now.
17-Mar-1999 MLR Removed /WAIT from caput, to be compatible with version 4.3 and later of the MCA record, which does not fire forward links until acquisition is complete.
28-Mar-1999 MLR Changed routine so it no longer pokes READ field. This assumes that someone else (typically a database) is periodically poking the READ field. The object initialization code now sets a monitor on the VAL field. Added New_flag output and CHECK_NEW keyword.

(See `epics_mca_define.pro`)




IDL EPICS Applications

- **mcaDisplay**
 - Full-featured program for displaying, controlling EPICS multi-channel analysers, including peak fitting
 - Uses epics_mca class library, and exports mca_display class, so it can be controlled by other IDL applications



IDL EPICS Applications

Data catcher and data viewer (Ben-Chin Cha)



Using EPICS from Visual Basic

- **ezca.dll** can be called directly from Visual Basic on Windows
- **ezca.bas** provides the interface

```
Public Const ezcaByte As Byte = 0
Public Const ezcaString As Byte = 1
Public Const ezcaShort As Byte = 2
Public Const ezcaLong As Byte = 3
Public Const ezcaFloat As Byte = 4
Public Const ezcaDouble As Byte = 5
```

```
Public Declare Function ezcaGet Lib "ezca.dll" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByRef data As Any) As Long
```




```
Public Declare Function ezcaPut Lib "ezca.dll" Alias "ezcaPutOldCa" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByRef data As Any) As Long
```


```
Public Declare Function ezcaPutString Lib "ezca.dll" Alias "ezcaPutOldCa" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByVal data As Any) As Long
```

```
Public Declare Function ezcaPutCallback Lib "ezca.dll" Alias "ezcaPut" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByRef data As Any) As Long
```

Using EPICS from Visual Basic

- Example: tomography data collection. VB used because it can easily control Roper's WinView program for the CCD detector





Python Applications for Beamline Control

Mark Rivers

Argonne National Laboratory



*A U.S. Department of Energy
Office of Science Laboratory
Operated by The University of Chicago*



Motivation

- **Replace IDL applications (e.g. MCA GUI) with Python so that other beamlines don't need to buy IDL**
- **Send users home with data and display/analysis programs that are free.**
 - They don't want to buy IDL.


Building Blocks

- **CaChannel from Geoff Savage for EPICS interface**
- **Tkinter and Pmw for GUIs**
- **Pmw.Blit for plots**
- **Numeric for arrays**

BltPlot: Enhancements to the Pmw.Blt.Graph widget

- **A standalone plotting widget, `BltPlot.BltPlot`. This widget has menus to:**
 - Configure all of the plot characteristics
 - Save and restore the plot settings and data
 - Print the plot to a file or printer
- **Methods (`BltPlot.plot` and `BltPlot.oplot`) to create a new plot, to overlay more data, etc.**
- **Designed to provide a rough emulation of the command line plotting capabilities of IDL.**

BltPlot



BltPlot: Enhancements to the Pmw.Blt.Graph widget

- **GUI routines to configure all of the plot characteristics, such as axes, markers, legends, etc.**
 - These routines work with any Pmw.Blt.Graph instance so they can be used from the standalone plotting widget in this package (BltPlot.BltPlot) or from any application that uses the Pmw.Blt.Graph widget
- **Routines to save and restore plot settings and data.**
- **Used in the mcaDisplay described later.**

BltPlot Dialogs

Configure Graph

File Settings

General Plot Configuration

Title:

X Label: Y Label:

Colors: Background: PlotArea: Crosshairs: Grid:

☐ Show Grid
☐ Show Legend

Trace Configuration

#	Label	Color	Thickness	Dashes	Symbol Type	Symbol Size
1	cos	<input type="text"/>	thin	<input type="text"/>	none	tiny
2	sin	<input type="text"/>	thin	<input type="text"/>	none	tiny

OK Apply Quit

BltConfigureAxis

X axis

Hide or show:

Title:

Title font: *-Helvetica-Medium-R-

Title color: #000000

Tick interval: 0.0

Minor ticks: 2

Tick length: 8

Tick font: *-Helvetica-Medium-R-

Axis color: #000000

Limits:

Minimum: 18.0

Maximum: 886.0

Scale:

OK Apply Quit

epicsPV: Subclasses Geoff Savage's CaChannel class

- If a PV name is given then the class constructor will do a `searchw()` by default.
- `setMonitor()` sets a generic callback routine for value change events. Subsequent `getw()`, `getValue()` or `array_get()` calls will return the value from the most recent callback, and hence do not result in any network activity or latency. This can greatly improve performance.
- `checkMonitor()` returns a flag to indicate if a callback has occurred since the last call to `checkMonitor()`, `getw()`, `getValue()` or `array_get()`. It can be used to increase efficiency in polling applications.

- **getControl()** reads the "control" and other information from an EPICS PV without having to use callbacks. In addition to the PV value, this will return the graphic, control and alarm limits, etc.
- **putWait()** calls **array_put_callback()** and waits for the callback to occur before it returns. This allows programs to use **array_put_callback()** synchronously and without user-written callbacks.

epicsMotor

- **Class library for EPICS motor record**
- **Methods:**
 - `move()`, `stop()`, `wait()`, `get_position()`, `set_position()`
- **Virtual attributes:**
 - `slew_speed`, `base_speed`, `high_limit`, `low_limit`, `done_moving`, `backlash`, `resolution`, etc.
- **Example use:**

```
from epicsMotor import *
m = epicsMotor('13LAB:m5')
m.move(10.)
m.wait()
m.get_position(dial=1, readback=1)
9.9609375
```


epicsScaler

- **Class library for EPICS scaler record**
- **Methods:**
 - start(), stop(), read(), wait(), get_label(), set_label()
- **Example use:**

```
>>> from epicsScaler import *
>>> s = epicsScaler('13LAB:scaler1')
>>> s.get_counts()
>>> s.read()
[0, 0, 0, 0, 0, 0, 0, 0]
>>> s.start(1.)
>>> s.wait()
>>> s.read()
[10000000, 0, 0, 0, 0, 0, 0, 0]
```

epicsLogger


GUI application for logging EPICS PVs to the screen and to a disk file



The screenshot shows the 'epicsLogger' application window. It has a menu bar with 'File' and 'Help'. Below the menu bar is a control area with 'Start' and 'Stop' buttons, an 'Update time (seconds)' field set to '40', and an 'Output file: None' label. The main area contains a table with 7 columns: 'Date and time', 'Load, Tons', 'Ram Ht, mm', 'room temperature, C', 'Sample Tc, mV', 'T (C)', and 'Volts'. The table has 16 rows of data, all from '26-Mar-2003' between '19:30:27' and '19:31:34'. The data values are as follows:

Date and time	Load, Tons	Ram Ht, mm	room temperature, C	Sample Tc, mV	T (C)	Volts
26-Mar-2003 19:30:27	0.43917	-25.16517	34.87248	-81.71421	-4939.15607	0.03511
26-Mar-2003 19:30:37	0.33953	-25.16598	34.87110	-1.42667	-87.51303	0.03511
26-Mar-2003 19:30:47	0.44808	-25.16619	34.86937	-31.10073	-2608.25114	0.03513
26-Mar-2003 19:30:57	0.29490	-25.16464	34.86736	-54.78670	-4882.51811	0.03514
26-Mar-2003 19:31:07	0.44804	-25.16483	34.86229	-61.13838	-5287.09796	0.03514
26-Mar-2003 19:31:17	0.29018	-25.16452	34.86673	-56.22646	-4987.64225	0.03513
26-Mar-2003 19:31:20	0.29018	-25.16452	34.86673	-56.22646	-4987.64225	0.03513
26-Mar-2003 19:31:22	0.29018	-25.16452	34.86673	-56.22646	-4987.64225	0.03513
26-Mar-2003 19:31:24	0.47972	-25.16535	34.86702	-53.39937	-4774.91081	0.03513
26-Mar-2003 19:31:26	0.47972	-25.16535	34.86702	-53.39937	-4774.91081	0.03513
26-Mar-2003 19:31:28	0.47972	-25.16535	34.86702	-53.39937	-4774.91081	0.03514
26-Mar-2003 19:31:30	0.47972	-25.16535	34.86702	-53.39937	-4774.91081	0.03514
26-Mar-2003 19:31:32	0.47972	-25.16535	34.86702	-53.39937	-4774.91081	0.03514
26-Mar-2003 19:31:34	0.28039	-25.16466	34.86693	-36.80112	-3194.25474	0.03514

mcaDisplay



mcaDisplay()

- Replacement for my IDL MCA display program
- Much nicer in many respects, since the Blt plot widget has many more east-to-use features than IDL's direct graphics
- Python object with callable methods, so it can be remotely controlled
- Device independent. It reads files and controls the “hardware_mca” class. “hardware_mca” can be subclassed for any hardware. Presently the EPICS MCA record is supported

Mca: Device-independent MCA class

- **Support classes: mcaROI, mcaCalibration, mcaElapsed, mcaPresets, mcaPeak, etc.**
- **Many methods: add_roi(), fit_background(), fit_peaks(), get_calibration(), set_calibration(), write_file(), read_file(), etc.**
- **Used as base class of epicsMca.**


epicsMca:

Subclass of hardwareMca, which is subclass of Mca

- All methods of Mca, plus start(), stop(), erase(), wait(), etc.
- Re-implements base class routines for set_calibration(), set_rois(), etc. to communicate with fields in the EPICS MCA record
- **Example use:**

```
from epicsMca import *
mca = epicsMca('13BMD:aim_adc1')
mca.erase()
mca.start()
mca.wait()
mca.write_file('test.001')
```

mcaPeakFit



44



Pioneering
Science and
Technology

Office of Science
U.S. Department
of Energy



mcaPeakFit Parameters

Peak Fit

File

Background Parameters

Exponent: 4 Top width: 0.0 Bottom width: 4.0 Tangent?: No Compression: 4

Fit background Re-plot background

Peak Fit Parameters

Initial energy calibration

Offset: 0.047064327 Slope: 0.017652117 Flag: Optimize

Initial FWHM calibration

Offset: 0.15 Slope: 0.0 Flag: Optimize

Peak parameters

	Label	Energy	Energy flag	FWHM	FWHM flag	Ampl. ratio
	Mn Ka	5.895	Fix	0.15	Global	0.0
Insert	Mn Ka,	5.8950,	0,	0.1500,	0,	0.0000
Append	Fe Ka,	6.4000,	0,	0.1500,	0,	0.0000
	Fe Kb,	7.0590,	0,	0.1500,	0,	0.0000
Delete	Ni Ka,	7.4720,	0,	0.1500,	0,	0.0000
	Cu Ka,	8.0410,	0,	0.1500,	0,	0.0000
	Zn Ka,	8.6310,	0,	0.1500,	0,	0.0000
Clear	Zn Kb,	9.5720,	0,	0.1500,	0,	0.0000
Sort	Ga Ka,	9.2430,	0,	0.1500,	0,	0.0000
Mark all						

Fit results file: [fit_results.txt](#)

Spreadsheet file: [fit_spreadsheet.txt](#)

Overwrite or append results and spreadsheet files: Append

Fit peaks Re-plot fit

mcaPeakFit Output

```
mcaPeakFit Results

*****
Fit of /usr/local/python/epics/basalt_xrf.002:3

Real time (seconds):      872.51
Live time (seconds):      872.51
Initial FWHM offset, slope:  0.150000  0.000000
Optimized FWHM offset, slope: 0.305848  0.002562
Initial energy offset, slope: 0.047064  0.017652
Optimized energy offset, slope: -0.046042  0.017947
# Iterations, function evals:      8      93
Chi squared:              5.92047e+09
Status code:              2
Time to fit:              3.246

      Peak      Energy      FWHM      Area      Background      Area/MDL      Area/Bkg
Mn Ka      5.895      0.3121      204198.9      19513.0      487.3      10.5
Fe Ka      6.400      0.3123      3079580.9      19196.0      7409.1      160.4
Fe Kb      7.059      0.3127      949900.2      17498.0      2393.7      54.3
Ni Ka      7.472      0.3129      0.0      15107.0      0.0      0.0
Cu Ka      8.041      0.3131      17554.6      15052.0      47.7      1.2
Zn Ka      8.631      0.3134      49345.0      15625.0      131.6      3.2
Zn Kb      9.572      0.3138      14772.8      16192.0      38.7      0.9
Ga Ka      9.243      0.3136      63076.2      15686.0      167.9      4.0

      OK
```

Mpfit

- **Generalized non-linear least squares data fitting**
- **Based on LMFIT from Minpack**
- **Originally translated to IDL by Craig Markwardt, I translated to Python**
- **Much faster and more accurate than the version provided in the Scientific Python package in `Scientific.Functions.LeastSquares`.**
- **Constraints, fixed parameters, analytic or numerical derivatives, etc.**
- **Used in `mcaPeakFit`**

Med: Device-independent multi-element detector class

- **Collection of Mca objects. Methods operate on all contained Mca objects. Example:**
 - `add_roi()`, `set_presets()`, `get_calibration()`, etc.
- **Used as base class of `epicsMed`.**

epicsMed

- **Subclass of Mca and Med**
- **All methods of Mca and Med, plus start(), stop(), erase(), wait(), etc.**
- **Re-implements base class routines for set_calibration(), set_rois(), etc. to communicate with fields in the EPICS MCA record**

- **Example use:**

```
from epicsMed import *
med = epicsMed('13GE1:med:', 16)
med.erase()
med.start()
med.wait()
med.write_file('test.001')
```

medDisplay

