

ALIAS Instrument Capabilities and Tracer Measurements

Dr. Lance E. Christensen
Dr. Christopher R. Webster

Atmospheric Laser Spectroscopy Group
Jet Propulsion Laboratory
September 19, 2005


ALIAS Characteristics

- Laser absorbance/2f detection, 80-m or 160-m pathlength.
- Four channels, broad spectral coverage:

1) Methane	3.4 μm	IC laser
2) Water Isotopes	6.7 μm	QC laser
3) Water Isotopes	6.7 μm	QC laser
4) Water	1.9 μm	Near-IR laser
- Regulated inlet flow - isokinetic.
- 16 L sec⁻¹ flow through instrument.
- Temperature of measured gas is regulated.


Inlet/airflow

Differential pressure measurement at inlet. Valve regulation of impedance.




Inlet/airflow

Isokinetic flow



Inlet/airflow


Orifice-plate Flowmeter confirms flow through instrument



Airflow/water signal.

No significant contamination from fore-optics.


Herriott cell: 80-160m. Near-IR laser: 0.14m



Airflow/water signal.


Wind-tunnel experiments, scale-model water experiments, and numerical analysis indicate:

- Virtually no end-to-end pressure difference.
- Gas is well mixed throughout the analysis cell.
- $1/e$ residence time is 1.7 sec.



Airflow/water signal.


Very little evidence of significant water outgassing.




Temperature control

All water particles are vaporized.


Spear: 500W heating




Temperaure control




Temperature is well regulated during flight.




Carbon Monoxide




CH₄



CH₄





July 3, 2005
descent






July 7, 2005
descent

Better Signal-to-Noise for water has been achieved.
Shielding and filtering have been improved.



ALIAS employs a novel semi-conductor laser



Dr. Rui Yang, (JPL)

Interband Cascade Lasers (ICL):

- Cover 2.5 to 5 μm spectral region.
- High power.
- Single-mode emission (DFB).
- High efficiency (5-20mA threshold).
- Progress towards room T devices.

Conclusions

- AVE-WIIF mission confirmed that ALIAS is:
 - Isokinetic
 - Isothermal
 - not contaminated by water outgassing
- Tracer sensitivities (3 sec):
 - CH₄ 50 pptv
 - N₂O 60 pptv
 - CO 300 pptv
- ALIAS has integrated and employed IC Lasers.